

## Polymerase Chain Reaction(PCR)

### Reagents

- 1) 5  $\mu$ l 10X Standard Taq Reaction buffer
- 2) 1  $\mu$ l 10 mM dNTPs
- 3) 1  $\mu$ l 10 mM Forward Primer
- 4) 1  $\mu$ l 10 mM Reverse Primer
- 5) Template DNA
- 6) 25  $\mu$ l Taq DNA Polymerase
- 7) Make up to 50  $\mu$ l Nuclease-free water

### Procedure

The protocol is composed of 3 basic steps: denaturation, annealing, and extension.

- 1) Preheat the thermocycler to 95°C.
- 2) Gently mix the reaction components. If necessary, centrifuge briefly to collect the liquid to the bottom of the tube.
- 3) If using a PCR machine without a heated lid, overlay the sample with mineral oil.
- 4) Transfer PCR tubes from ice to the preheated PCR machine (95°C) for the initial denaturation step.
- 5) Begin thermocycling:
  - a) Initial denaturation: 95 °C, 30 seconds
  - b) 30 cycles of: 95 °C, 15-30 seconds; 45 °C – 68 °C, 15 – 60 seconds; 68 °C, 1 minute/kb
  - c) Final extension: 68 °C, 5 minutes
  - d) Hold: 4 – 10 °C
- 6) Validate the final reaction mixture by running a fraction on an agarose gel with an appropriate molecular marker. Directly sequence the amplified product or carry out restriction enzyme digests.

### Recommendations

- 1) MgCl<sub>2</sub>: For most Taq DNA Polymerase a MgCl<sub>2</sub> concentration of 1.5 – 2.0 mM is optimal. The final concentration in standard 1X Standard Taq Reaction Buffer is 1.5 mM. If

another concentration is optimal, a reaction buffer that does not contain  $MgCl_2$  can be used and the  $MgCl_2$  can be added separately.

- 2) Additives: DMSO and formamide can be added to improve difficult targets such as GC-rich sequences.
- 3) Template DNA: Use high quality, purified DNA. For a 50  $\mu$ l reaction, recommended amounts are:
  - a) Genomic DNA 1 ng – 1  $\mu$ g
  - b) Plasmid or viral DNA 1 pg – 1 ng
- 4) Taq DNA Polymerase: 1.25 units/50  $\mu$ l reaction is recommended, however this can vary from 0.25 – 2.5 units/50  $\mu$ l reaction.
- 5) Denaturation: A longer initial denaturation step (2 – 4 minutes at 95 °C) can be employed for difficult templates such as GC-rich sequences. For colony PCR 5 minutes at 95 °C is recommended. During thermocycling 15 – 30 seconds at 95 °C is recommended.
- 6) Annealing: The temperature is dependent on the  $T_M$  of the primer pair (usually 45 °C – 68 °C). Optimise by carrying out temperature gradient PCR starting 5 °C below the  $T_M$ .
- 7) Extension: The recommended temperature is 68 °C with the time generally 1 minute per kb. A final extension of 5 minutes at 68 °C is recommended.
- 8) Cycle number: Generally, 25 – 35 cycles yields sufficient product. Up to 45 cycles can be employed to detect low copy number targets.