

## MEK2 Rabbit Polyclonal Antibody

### Catalog #: EAB11894

| Host/Isotype | Clonality  | Applications  | MW (kDa) | Reactivity        |
|--------------|------------|---------------|----------|-------------------|
| Rabbit IgG   | Polyclonal | WB, IHC-P, IF | 44       | Human, Mouse, Rat |

### Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|                                              |             |
|----------------------------------------------|-------------|
| <b>WB</b> (Western Blotting)                 | 1:1000-5000 |
| <b>IHC-P</b> (Immunohistochemistry-Paraffin) | 1:100-500   |
| <b>IF</b> (Immunofluorescence)               | 1:50-300    |

### Product Information

|                       |                                                                                                                                           |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Conjugate</b>      | Unconjugate                                                                                                                               |
| <b>Specificity</b>    | MEK2 Rabbit Polyclonal Antibody detects endogenous levels of MEK2 protein.                                                                |
| <b>Purification</b>   | Affinity purification                                                                                                                     |
| <b>Concentration</b>  | 1mg/ml                                                                                                                                    |
| <b>Format</b>         | Liquid                                                                                                                                    |
| <b>Formulation</b>    | In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol                                                                  |
| <b>Shipping</b>       | Gel Pack                                                                                                                                  |
| <b>Storage</b>        | Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles.<br>Aliquots may be stored at +4°C for 1-2 weeks |
| <b>UniProt ID</b>     | <a href="#">P36507</a>                                                                                                                    |
| <b>Entrez-Gene Id</b> | <a href="#">5605</a>                                                                                                                      |

### Product Description

The protein encoded by this gene is a dual specificity protein kinase that belongs to the MAP kinase kinase family. This kinase is known to play a critical role in mitogen growth factor signal transduction. It phosphorylates and thus activates MAPK1/ERK2 and MAPK2/ERK3. The activation of this kinase itself is dependent on the Ser/Thr phosphorylation by MAP kinase kinase kinases. Mutations in this gene cause cardiofaciocutaneous syndrome (CFC syndrome), a disease characterized by heart defects, cognitive disability, and distinctive facial features similar to those found in Noonan syndrome. The inhibition or degradation of this kinase is also found to be involved in the pathogenesis of Yersinia and anthrax. A pseudogene, which is located on chromosome 7, has been identified for this gene.

**For Reserch Use Only. Not For Use In Diagnostic Procedures**