

PPP6C Rabbit Polyclonal Antibody

Catalog #: EAB13104

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB	35	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB(Western Blotting) 1:500-2000

Product Information

Conjugate	Unconjugate
Specificity	PPP6C Rabbit Polyclonal Antibody detects endogenous levels of PPP6C protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol.
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks.
UniProt ID	Q00743
Entrez-Gene ID	5537

Product Description

PPP6C also known as PP6C, PP6 (protein phosphatase 6) is a 305 amino acid cytoplasmic protein that belongs to the PPP phosphatase family. Ubiquitously expressed, PP6 is component of a signaling pathway regulating cell cycle progression in response to IL-2 receptor stimulation and is involved in suppressing inflammatory responses by specifically down-regulating TR4.

For Reserch Use Only. Not For Use In Diagnostic Procedures