

## Phospho-ALK (Tyr1278) Rabbit Polyclonal Antibody

### Catalog #: EAB13210

| Host/Isotype | Clonality  | Applications | MW (kDa) | Reactivity        |
|--------------|------------|--------------|----------|-------------------|
| Rabbit IgG   | Polyclonal | WB           | 176      | Human, Mouse, Rat |

### Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

**WB**(Western Blotting) 1:500-2000

### Product Information

|                       |                                                                                                                                        |
|-----------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Conjugate</b>      | Unconjugate                                                                                                                            |
| <b>Specificity</b>    | Phospho-ALK (Tyr1278) Rabbit Polyclonal Antibody detects endogenous levels of ALK only when phosphorylated at Tyr1278.                 |
| <b>Purification</b>   | Affinity purification                                                                                                                  |
| <b>Concentration</b>  | 1mg/ml                                                                                                                                 |
| <b>Format</b>         | Liquid                                                                                                                                 |
| <b>Formulation</b>    | In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol                                                               |
| <b>Shipping</b>       | Gel Pack                                                                                                                               |
| <b>Storage</b>        | Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks |
| <b>UniProt ID</b>     | <a href="#">Q9UM73</a>                                                                                                                 |
| <b>Entrez-Gene Id</b> | <a href="#">238</a>                                                                                                                    |

### Product Description

This gene encodes a receptor tyrosine kinase, which belongs to the insulin receptor superfamily. This protein comprises an extracellular domain, an hydrophobic stretch corresponding to a single pass transmembrane region, and an intracellular kinase domain. It plays an important role in the development of the brain and exerts its effects on specific neurons in the nervous system. This gene has been found to be rearranged, mutated, or amplified in a series of tumours including anaplastic large cell lymphomas, neuroblastoma, and non-small cell lung cancer. The chromosomal rearrangements are the most common genetic alterations in this gene, which result in creation of multiple fusion genes in tumourigenesis, including ALK (chromosome 2)/EML4 (chromosome 2), ALK/RANBP2 (chromosome 2), ALK/ATIC (chromosome 2), ALK/TFG (chromosome 3), ALK/NPM1 (chromosome 5), ALK/SQSTM1 (chromosome 5), ALK/KIF5B (chromosome 10), ALK/CLTC (chromosome 17), ALK/TPM4 (chromosome 19), and ALK/MSN (chromosome X).

**For Reserch Use Only. Not For Use In Diagnostic Procedures**