

Phospho-FRS2 (Tyr196) Rabbit Polyclonal Antibody

Catalog #: EAB13297

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB	65	Human, Mouse

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB(Western Blotting) 1:500-2000

Product Information

Conjugate	Unconjugate
Specificity	Phospho-FRS2 (Tyr196) Rabbit Polyclonal Antibody detects endogenous levels of FRS2 only when phosphorylated at Tyr196.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	Q8WU20
Entrez-Gene Id	10818

Product Description

FRS2 (also designated SNT or p90) is a lipid-anchored docking protein that becomes tyrosine phosphorylated in response to FGF or NGF stimulation and subsequently binds to GRB2/Sos complexes. The GRB2 adapter protein links receptor tyrosine kinases to the Ras/MAPK signaling pathway but does not interact directly with FGF receptors. FRS2 thus provides a link between activation of FGF and NGF receptors and the Ras/MAPK pathway. FRS2 contains four Grb2 binding sites, a myristylation sequence and a PTP domain. Myristylation of FRS2 is essential for membrane localization, tyrosine phosphorylation, GRB2/Sos recruitment and MAPK activation. The function of FRS2 in FGF receptor signaling is analogous to that of IRS1 in response to insulin receptor stimulation.

For Reserch Use Only. Not For Use In Diagnostic Procedures