

Rrn3 Rabbit Polyclonal Antibody

Catalog #: EAB14736

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, IHC-P, ELISA	70	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB (Western Blotting)	1:500-2000
IHC-P (Immunohistochemistry-Paraffin)	1:50-300
ELISA (Enzyme-linked Immunosorbent Assay)	1:5000-20000

Product Information

Conjugate	Unconjugate
Specificity	Rrn3 Rabbit Polyclonal Antibody detects endogenous levels of Rrn3 protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol.
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks.
UniProt ID	Q9NYV6
Entrez-Gene ID	54700

Product Description

Enables RNA polymerase I core promoter sequence-specific DNA binding activity. Predicted to be involved in transcription initiation from RNA polymerase I promoter. Predicted to act upstream of or within several processes, including DNA-templated transcription, initiation; negative regulation of intrinsic apoptotic signaling pathway by p53 class mediator; and regulation of transcription, DNA-templated. Located in nucleolus.

For Reserch Use Only. Not For Use In Diagnostic Procedures