

γ-Catenin Rabbit Polyclonal Antibody

Catalog #: EAB21228

Host/Isotype	Clonality	Applications	MW (kDa)	Reactivity
Rabbit IgG	Polyclonal	WB, IHC-P, IF/ICC	82	Human, Mouse, Rat

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

WB (Western Blotting)	1:5000-20000
IHC-P (Immunohistochemistry-Paraffin)	1:500-2000
IF/ICC (Immunofluorescence/Immunocytochemistry)	1:100-500

Product Information

Conjugate	Unconjugate
Specificity	γ-Catenin Rabbit Polyclonal Antibody detects endogenous levels of γ-Catenin protein.
Purification	Affinity purification
Concentration	1mg/ml
Format	Liquid
Formulation	In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol
Shipping	Gel Pack
Storage	Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks
UniProt ID	P14923
Entrez-Gene ID	3728

Product Description

The catenins, α, β and γ, are proteins which bind to the highly conserved, intracellular cytoplasmic tail of E-cadherin. Together, the catenin/cadherin complexes play an important role mediating cellular adhesion. γ-catenin, also known as plakoglobin, binds with α-catenin and N-cadherin. It has been shown that the transmembrane phosphatase PTPμ associates with catenin/cadherin complexes and may regulate complex signaling.

For Reserch Use Only. Not For Use In Diagnostic Procedures