

MMP-12 Rabbit Monoclonal Antibody

Catalog #: EAB21360

| Host/Isotype | Clonality | Applications | MW (kDa) | Reactivity |
|--------------|------------|---------------------------|----------|-------------------|
| Rabbit IgG | Monoclonal | WB, IP, IHC-P, IF/ICC, FC | 54 | Human, Mouse, Rat |

Applications Dilutions

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|--|------------|
| WB (Western Blotting) | 1:500-2000 |
| IP (Immunoprecipitation) | 1:10-100 |
| IHC-P (Immunohistochemistry-Paraffin) | 1:50-200 |
| IF/ICC (Immunofluorescence/Immunocytochemistry) | 1:50-200 |
| FC (Flow Cytometry) | 1:10-100 |

Product Information

| | |
|-----------------------|---|
| Conjugate | Unconjugate |
| Specificity | MMP-12 Rabbit Monoclonal Antibody detects endogenous levels of MMP-12 protein. |
| Purification | Affinity purification |
| Concentration | 1mg/ml |
| Format | Liquid |
| Formulation | In PBS, pH 7.4, Containing 0.02% sodium azide, 0.5% BSA and 50% Glycerol |
| Shipping | Gel Pack |
| Storage | Store at -20°C least 1 year from the date of shipment. Avoid repeated freeze/thaw cycles. Aliquots may be stored at +4°C for 1-2 weeks |
| UniProt ID | P39900 |
| Entrez-Gene Id | 4321 |

Product Description

This gene encodes a member of the peptidase M10 family of matrix metalloproteinases (MMPs). Proteins in this family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. The encoded preproprotein is proteolytically processed to generate the mature protease. This protease degrades soluble and insoluble elastin. This gene may play a role in aneurysm formation and mutations in this gene are associated with lung function and chronic obstructive pulmonary disease (COPD). This gene is part of a cluster of MMP genes on chromosome 11.

For Reserch Use Only. Not For Use In Diagnostic Procedures